

# Microfabricated Organic Analyzer (MOA) for *in situ* Exploration of Mars and other Solar Bodies

- Prof. Richard A. Mathies, Chemistry Department, UC Berkeley, microfabricated devices and biochemical analysis
- Prof. Jeff Bada, Scripps Institution of Oceanography, UCSD, amino acid analysis and astrobiology
- Dr. Frank Grunthaner, Jet Propulsion Laboratory, instrument design, operation and flight engineering

A wide-angle photograph of the Martian surface, showing a vast, flat, reddish-brown landscape under a hazy orange sky. In the lower-left foreground, a small portion of a rover is visible, including a white cylindrical component and some mechanical parts. The terrain is covered in fine-grained soil and small, dark rocks.


# Project Goals

- Complete brass-board development of microfabricated capillary electrophoresis (CE) chip and instrument for amino acid analysis
- Integrate the microchip CE system with MOD sampling system to form the Mars Organic Analyzer (MOA)
- Perform field tests of the MOA in Mojave and three Mars-like Atacama sites
- Document maturation, integration and field operation of two MIDP, PIDDP and ASTID-derived instruments
- Enhance TRL of MOA through Mojave and Atacama field testing
- Increase our understanding of limits and constraints of life in extreme environments
- Critically define identity of and sensitivity requirements of potential biomarkers.

# Amino Acid Composition and Chirality Analysis

## Potential bioorganic signatures:


- Large biomolecules likely degraded by oxidizing surface environment of Mars
- Amino acids have a longer lifetimes in dry, harsh conditions
- Amino acids have been found in meteorites
- Amino acid chirality is indicative of origin:
  - Racemic mixture – abiotic origin
  - Non-racemic mixture – biological origin



L-alanine

## Why *in situ* analysis:

- Significant contamination of Meteorites found on earth by terrestrial sources of life
- Sample return missions are more technologically challenging, costly and time consuming




Allan Hills Meteorite



# Mars Organic Analyzer (MOA) Concept

- Soil samples collected and deposited into Mars Organic Detector (MOD)
- MOD sublimates amino acids onto cold finger coated with fluorescamine<sup>1</sup>
- Fluorescamine-labeled amino acids analyzed for composition and chirality via microchip Capillary Electrophoresis



<sup>1</sup> Kminek, et al., *Planetary & Space Science* 2000, 48, 1087-1091.



# Composition and Chirality Analysis by CE



- Electroosmotic flow (EOF) sweeps all molecules to the cathode

- CE separates amino acids based on charge/size ratio giving composition information

- Cyclodextrins included in running buffer provide enantiomeric resolution of amino acids




$$K_L \approx K_D$$




L-alanine

D-alanine




g-cyclodextrin

# Established Separation and Lab-Based Detection System



## Analysis of Fluorescein-labeled Amino Acids



## Murchison Meteorite Glu and Asp D/L values


Amino Acid	HPLC	Microchip CE
<b>Glu</b>		
Exterior	0.3 - 0.1	0.33 ± 0.04
Interior	0.7 - 0.1	0.65 ± 0.07
<b>Asp</b>		
Exterior	0.3 - 0.1	0.21 ± 0.03
Interior	0.3 - 0.1	0.30 ± 0.06

Hutt, L. D.; Glavin, D. P.; Bada, J. L.; Mathies, R. A. *Analytical Chemistry* 1999, 71, 4000-4006.


# Separation of Fluorescamine-labeled Amino Acids

## Advantages of Fluorescamine:

- Fluorogenic reagent
- Reaction time ~1 min
- ~ 50 nM LOD attainable
- Reagent used in MOD




## Separation of Mars 7 Standard labeled with Fluorescamine



A. M. Skelley and R. A. Mathies, *J. Chromatogr. A* 2003, 1021, 191-199.




# Atacama Desert as Martian Analog Site



- Chilean Atacama Desert is one of the driest sites on the planet ( $<0.5$  mm  $H_2O$ /year).
- The transect from Lat  $24^\circ$  to  $28^\circ$  South at  $69.5^\circ$  West has been extensively studied.
- Some areas have unusual surface oxidation chemistry and organic soil concentrations at lab blank levels. Other areas show readily detected microbial and higher life forms.

# Analysis of Atacama Soil Extracts


## Composition



### Blank-corrected Concentrations (weight/weight):

Val =  $0.034 \pm 0.009$  ppm      Ala/Ser =  $0.32 \pm 0.07$  ppm  
 Gly =  $0.18 \pm 0.03$  ppm      Glu =  $0.14 \pm 0.02$  ppm  
 Asp =  $0.094 \pm 0.004$  ppm

## Chirality





### Blank-corrected Enantiomeric Ratios:



Mars 7 Standard:      Atacama Sample # 20  
 D/L Glu =  $1.10 \pm 0.02$       D/L Glu =  $0.22 \pm 0.02$   
 D/L Asp =  $0.97 \pm 0.02$       D/L Asp =  $0.16 \pm 0.02$



# Portable Microchip CE System - Schematic




**W.H. Grover, A.M. Skelley, C.N. Liu, E.T. Lagally, R.A. Mathies, *Sens. Actuators B* 89 (2003) 325.**






# Portable Microchip CE Instrument




# Comparison of Lab and MOA Systems

**Concentration of each AA = 2.2 mM**






**Concentration of each AA = 444 nM**



**• MOA system shows superior sensitivity and comparable separation efficiency**

# MOD + CE = MOA







# Summary

- **Amino acid composition and chirality is an ideal means for organic biomarker detection on Mars.**
- **Microfabricated CE instrument provides a demonstrated means for sensitive amino acid composition and chirality analysis.**
- **Portable CE instrument has identical separation efficiency and greater sensitivity than standard lab systems.**
- **Integration of CE with MOD to make MOA will provide sensitive analysis of amino acids in Martian soil.**
- **Field tests in Mojave and Atacama Deserts are planned as a critical test of technology readiness and analysis capabilities.**
- **Microchip is a powerful platform for preparation other analytes from other sources for many types of analyses.**


# Mars Organic Laboratory (MOL): Beyond amino acids



**MOD**




**Sub-Critical Extractor**



**Cryobot**

## MOL Schematic - ASTID Project



1. A.M. Skelley et al., *SPIE: Proceedings of the In-Situ Instrument Technologies Meeting*, 4878 (2002) 59.

2. E.T. Lagally et al., *Lab-on-a-Chip*, 1 (2001) 107.

3. E.T. Lagally et al., *Anal. Chem.*, submitted (2004).

Integrated PCR with CE analysis of nucleic acids

Portable bacterial detection and typing instrument

# Acknowledgements

- Alison Skelley, Jim Scherer and Will Grover
- ASTEP - *Microfabricated Organic Analyzer (MOA) for in situ Exploration of Mars and Other Solar Bodies* PI's Mathies, Grunthaner and Bada, 1/04 -1/07.
- ASTID - *Integration of a Micro-chip Amino Acid Chirality Detector into the Mars Organic Detector*, subcontract to UCB from UCSD, PI Bada, 5/02 - 5/04.
- JPL Instrument Development Funds
- NSERC Fellowship to Alison Skelley
- UC Berkeley Microfabrication Laboratory