



An Integrated/Modular Sample Handling and Biotechnology Life Detection Instrument for Solar System Exploration

PI – Andrew Steele

Geophysical Laboratory
Carnegie Institution of Washington

Name	Affiliation	Status in Proposal
Andrew Steele	CIW	PI
Jan Toporski	CIW	Co-I
Marilyn Fogel	CIW	Co-I
Wesley Huntress	CIW	Co-I
Jake Maule	CIW/JSC	Co-I
Victor Parro Garcia	CAB	Co-I
Carlos Briones	CAB	Co-I
Felipe Gómez	CAB	Collaborator in Spain
David Fernández-Remolar	CAB	Collaborator in Spain
Olga Prieto	CAB	Collaborator in Spain
Javier Gómez	CAB	Collaborator in Spain
Juan Pérez-Mercader	CAB	Collaborator in Spain
Fernando Rull	CAB	Collaborator in Spain
Norm Wainwright	Marine Biology Lab	Co-I
Deborah Kelley	Washington University	Co-I
John Baross	Washington University	Co-I
David McKay	NASA JSC	Co-I
John Lindsay	NASA JSC	Co-I
Jackie Allen	NASA JSC	Collaborator
Kennda Lynch	LMSO NASA JSC	Collaborator
Mary Schweitzer	Montana State University	Co-I
James Schmitt	Montana State University	Co-I
Judson Hedgecock	Oceaneering Space Systems	Co-I
Yosef Bar-Cohen	JPL's NDEAA Lab	Co-I
Seth Pincus	New Orleans Children's Hospital	Co-I
Huw Rowlands	Blue Sky Project Management	Consultant
Pascal Lee	SETI Institute	Co-I
Charlie Cockell	SETI Institute	Co-I
Wendell Mendell	NASA JSC	Collaborator
M. Wuttke	Department of Natural History, Mainz	Collaborator
Ken Nealson	USC	Collaborator

Philosophy

To define a non-rover mobile platform containing remote sample preparation capabilities including; sample acquisition, drilling or sub coring and secondary sample processing that can be used by any instrument for testing in any Mars relevant environments.

To make this system light, flexible, easy to integrate to and available to any instrument team for use in their testing.

To show this by using existing instrumentation combined with instruments under development in several environments.

To test in field and laboratory off the shelf instrumentation to measure viable and fossil organic biosignatures in all environments to verify instruments under development. Includes detection sensitivity, performance in the field and contamination issues.

Science Objectives

S1	Define suitable probes (e.g., DNA, aptamers, antibodies) for several terrestrial field sites through assay in the laboratory of returned samples from these sites.
S2	Identify optimal techniques for <i>in-situ</i> extraction of important biomarkers located in rocks, soil, and liquid samples.
S3	Test DNA/protein micro-array technology, PCR, fluorescence, and enzyme based life detection instruments in the above environments.
S4	Integrate these with the life detection instruments including SOLID (DNA/protein chip technology), PCR, LAL enzyme assay, ATP fluorescence, microscopy and fluorometric analysis.
S5	Test regimes to ensure correct mineral and pH buffering for each technique mentioned above.
S6	Study <i>in-situ</i> terrestrial microbial populations associated with solar system relevant environments
S7	Define the detection sensitivity of a suite of biotechnology instruments in the field.

Technology Objectives

T1	Develop an integrated robotic biotechnology platform to acquire samples and detect life in Astrobiology relevant environments. Use of Ultrasonic corer to process drill cores.
T2	Develop requirements and a design for an autonomous sample acquisition and handling system capable of obtaining and processing soil, rock, and liquid samples from a variety of relevant environments.
T3	Integrate custom and off-the-shelf biotechnology instruments with a custom microfluidics system to produce a compact suite of instruments.
T4	Refine the instrument suite design based upon the results of field tests.

Field Campaign Objectives

F1	Integrate laboratory testing with the first phase of field testing in a “safe”, “close-to-services” environment.
F2	Determine the detection sensitivity of a suite of biotechnology instruments in relevant environments. Coordinate logistics including shipping expensive technology, testing and assembling a field lab, and using video cons and interactive web page design for debriefing and dissemination of “lessons learned.”
F3	Test Instrumental capabilities in the field including a) Communications technology; b) Robotic aspects; c) Power requirements and capabilities; d) Base plate flexibility and functionality; e) Computer integration; and f) System integration in the field; g) Test biotechnology instruments.
F4	Test Environmental effects in the field to determine whether the instrument can withstand and function up to specifications in different conditions of a) Landscape form; b) Substrate form (soil vs. rock); c) Ambient Temperature; d) Relative Humidity (Aqueous to Dry); e) Salinity; and f) Weather, e.g., Wind, Dust, or Precipitation.

INITIAL PROPOSED FIELD SITES

- v Home Station, NASA JSC
- v Rio Tinto, Spain
- v Hot Springs (Yellowstone Caldera)
- v Deep Sea Hydrothermal Systems (without described platform)
- v Enspel Formation, Enspel, Westerwald, Germany (Fossil Life Relevant)
- v Australian Desert Soils and Precambrian Rocks (Marble Bar and the Apex Chert)
- v Haughton Crater, Devon Island, Canada

Each field location had a different PI and organizing scientists and was designed to test different Mars relevant scenarios over a three year period.

Technology







Table 2	Mass	Dimensions	Power (W)
SimRP Platform	41.4kg	104.2cm x 34.3cm x 94.0cm	Nil
Robotic Arm	7.5 Kg	Reach 101 cm	25-75
Drill	> 1 Kg	25cm inches drill depth of 12-13 cm	5-7
	Voltage (VDC)	Current (mA)	Power (W)
Panasonic Camera	12	780	9.36
NI Compact FieldPoint System	20	800	16



Platform




Where possible the platform was designed to use off the shelf components










- Direct connectivity to sensors and actuators
- 8 and 16-channel modules; individually configurable channels
- Hot-swappable and autoconfigurable
- Programmable power-up states
- -40 to 70 °C operating range, 50G shock 5G vibration
- 8 Analog and Digital I/O Modules for Compact FieldPoint
- LabVIEW Real-Time embedded controllers for control, measurement, and signal processing
- Operate as stand-alone embedded real-time controllers or PC-based distributed I/O Ethernet interface (connection via satellite phone).
- Embedded Web and file servers with remote panel user interface with real time updates
- Nonvolatile memory for embedded program storage and data logging
- Onboard DRAM memory for embedded program operation
- 2 RS-232 and 1 RS-485 serial ports for connection to peripherals

System Management

Date and Time

Main Power Bus Voltage



ABLISSE Control and Data Handling

Device Management

Robot Arm

US Drill

Microscope/Camera

Microfluidics

ATP Luminometry

LAL

PCR

SPISE3

Test Parameters

Method

Endpoint-Kinetic

Endpoint

Kinetic-X

OperatorID

Cartridge

Data

Start Test!



Download File Name

%

Coefficients

Coeff1

0.00

Coeff2

0.00

Coeff3

0.00

Coeff4

0.00

EU Min/Max

0.00


EU Max

0.00

ECS Slope

0.00

LAL Data



Iteration

0

SHUTDOWN
PROGRAM

ABLISSE Control and Data Handling

System Management

Date and Time

5/1/2003 12:27:18

Main Power Bus Voltage



Comm Link

Transmit in Progress

Device Management

Robot Arm

US Drill

Microscope/Camera

Microfluidics

ATP Luminometry


LAL

PCR

SPISE3

Health and Status

Drill Temperature (C)



Overheat
Warning



Drill Head Pressure (psi)



Overpressure
Warning



Drill Depth (mm)

Drill Rate (mm/s)



Time (s)

Calculated Drill
Rate (mm/s)

0.00

Drill Limit



Command

Drill Motor

On



Off

Drill Depth Limit (mm)



0.00

Drill SHUTDOWN!



Drill SHUTDOWN Override



Iteration

344



SHUTDOWN
PROGRAM

Bulk Sample Handling and Delivery System

Provide Liquid sample handling.

Initially provided by CAB – SOLID instrument or MASSE instrument
Neither instrument are currently ready for field deployment and therefore a separate sample handling has been designed and is currently being fabricated.


- Wetting and mixing of the initial soil sample
 - pH Buffering and Metered Reagent Delivery
 - Heated Incubation (if needed)
 - Macro-filtration to remove particulates
 - Micro-filtration and pre-concentration of sample.
 - Metered microfluidic delivery of sample to each analysis instrument
-
- Initially until completed a manual fluid interface would be used.



INSTRUMENTS

The instrument suite can be used both as a primary biological detection suite or to verify the measurements made on other instruments

Top arrow shows Optical Chamber
where the reaction is measured



Bottom arrow shows 25 μ L Sample Chamber



Limulus Amebocyte Lysate
(LAL)

Charles River Endosafe unit, N
Wainwright MBL.

Detects Endotoxin (LPS) and is
being used in Planetary Protection
projects including having
swabbed the MER rovers to detect
bio-load of Gram Negative
bacteria. Being developed for
space flight.


Detection Sensitivity – 1 Cell (if
you can get it into the machine)
Extraction



ATP luminometer

Reads amount of ATP within seconds. Records metabolic activity in solid, liquid and swabbed samples. Being developed for space flight.

Detection sensitivity $10^2 - 10^3$ cells.





Mobile PCR

Complete mobile molecular laboratory. Enables extraction, PCR, Gel loading, Gel imaging within a single case.

PCR is as sensitive as extraction method.

Protocols have been optimized to ensure no contamination whilst in the field.

Several extraction methods have been studied including MoBio, ($10^4 - 10^6$ cells), DNA pure (10^6) and FTA paper (10^2 cells in solution).



DG2 Digital Microscope Scalar scopes.

Lenses contain light source and are made to be in focus upon contact with the surface of interest.

Includes, x30 – x 1000 lenses. UV, bore hole, endoscope attachments.

Others in the pipe line or in the proposal-

SPISE3 – Pan Conrad.

Hand held Capillary Electrophoresis

Hand Held PCR instruments

RT-PCR Instruments

Portable Raman spectrometer

Deployment

LAL – Planetary protection, Chesapeake bay, CIW, Svalbard (AMASE), rio Tinto.

ATP – Belize, Enspel, Chesapeake bay, Svalbard, rioTinto.

PCR – Chesapeake bay, Svalbard, rio Tinto, MSFC

Microscopy - Belize, Enspel, Chesapeake bay, Svalbard, rioTinto, Taiwan.

To come – supporting P Conrad ASTEP, Taiwan (deep drilling), Italy (drilling volcanoes), Svalbard.

Chesapeake Bay


Test contamination and field trip protocols in a near service environment.



A

PCR – 16S rDNA, Archaea specific, Ammonia Oxidase and Nitrate reductase.

All negatives were negative, including extraction blanks.






B


LAL and ATP worked well, results showed $10^6 - 10^8$ cells within the ore taken for PCR.

Arctic Mars Analogue Svalbard Expedition (AMASE)



Why Svalbard ?






Science Goals

Life detection and characterization in Mars analogue environment

How did the Carbonate globules form – biotic V abiotic

How do the organisms in the warm springs contribute to the formation of the carbonate terraces

Troll Springs




Primer Pair	Target
B27F – 1492R	Generic Eubacteria
PufM.557F – PufM.750R	Specific – purple sulfur bacteria*
GS.619F – GS1144R	Specific – Green sulphur bacteria 16S rDNA
CFX.856F – CFX.1240R	Specific - Green Non sulphur bacteria 16S rDNA
HB.418F – HB.1159R	Specific - Heliobacteria 16SrDNA
NirK.F – NirK.5R	Specific – Nitrate reductase
Amoa.1F – Amoa.2R	Specific – Ammonia Oxidase
APS1.Fb – APS1.Rb	Specific - APS reductase
FDH.F – FDH.R	Specific – Formate Dehydrogenase
McrA.1F – McrA.1R	Specific – Methyl coenzyme M reductase (methanogens)
Nu-SSU-0817 – Nu-SSU-1536	Generic - Eukaryotes

Pri mer / Area	Eub	Euk	Pb	GS	GN S	HB	Nr	Ao	APS	FD H	Mcr A
1	++ +	-	+	+	-	++	+	-	++	++	-
2	++ +	-	++ +	++ +	+	-	+	-	++ +	+	++ +
3	++ +	-	++ +	-	-	-	-	-	++	+	-
4	++ +	-	++ +	-	-	-	-	-	-	+	-
5	++ +	-	++ +	-	-	-	-	-	-	-	-
7	++ +	+	++ +	-	-	-	-	-	++ +	+	-
Blank	+	-	-	-	-	-	-	-	-	-	-
PCR blank	-	-	-	-	-	-	-	-	-	-	-

Eub – Eubacterial; Euk – Eukaryotic; Pb – Purple Bacteria; GS – Green sulphur bacteria; GNS – Green Non sulphur bacteria; HB – Heliobacteria; Nr – Nitrate reductase; Ao – Amoonia oxidase; APS – APS reductase; FDH – Formate dehydrogenase; McrA – Methanogenesis. - - negative; + - weakly positive, ++ - medium positive (similar band fluorescence to ladder) +++ - strongly positive.

In – Field ATP and LAL




Area	1	2	3	4	5	6	7
ATP	6140	229	5644	2383	NM	1730	NM
LAL	>0.155	<0.100	0.340	0.648	<0.100	NM	0.463


Lab verification of these results (and PCR) is underway. Metabolism and cell mass do correlate in certain areas but do not in others, however all detect life.

Heater on LAL assay needed and has been improved


Carbonate terraces – Troll springs




ATP – Analysis of cryptoendolithic communities

















- Sample 1 – Free floating filaments emerging with the spring water.
 - Sample 2 – Biofilm filaments along the wall of the spring
 - Sample 2b - Biofilm filaments along the wall of the spring
 - Sample 3 – Rubbing of the FTA paper along green (cyanobacterial biofilm) on the wall of the cave above the spring.
- Blank – Single disk of unused FTA paper that underwent same extraction protocol as the samples outlined above.




	Eub	Euk	APS	Pb	GS	McrA
1	+++	-	++	-	-	-
2	+++	-	+	-	-	+
2b	+++	-	+	-	-	-
3	-	-	-	-	-	-
Blank	-	-	-	-	-	-

Carbonates









ASTEP Conclusions

We have defined an instrument suite that works in labs and in field for the characterization of microbes

Addressed contamination issues

Defined a simple platform for instrument testing

Looking to collaborate with any and all who wish to use such an instrument suite for verification of their instruments or as a stand alone package

Future

Other instruments and deploy array based instruments when possible

Check other life detection strategies

Refine platform design

Svalbard -----

Done 3 environments, several more this year with Pan

Acknowledgements

All CIW ASTEP 03 / 04 collaborators and Co_Is

Pan Conrad and SPISE 3 team

Wes Huntress, James Hall

Maia Schweizer, Hans Amundsen and the AMASE
team

MASSE team

MARTE team